



MED-LOGICS, INC.

Advanced Options for M2 Users

Accurate
Durable
Epi-friendly



OEM



Loose fit between blade assembly and head results in an unpredictable flap thickness.

MED-LOGICS



Tighter tolerance resulting in more accurate flap thickness.

OEM









Accuracy between heads can deviate and produce inaccurate flap thicknesses.

MED-LOGICS



Every ML7050 M Head holds ± 5 microns in the critical dimension insuring greater flap thickness accuracy.

CustomFlap™ Options

CustomFlap™ Thickness Options (in microns)	 MINUS 30 CLIP	 MINUS 20 CLIP	 MINUS 10 CLIP	 PLANO CLIP	 PLUS 10 CLIP	 PLUS 20 CLIP
ML7050 R100M Head (100 Micron)	70	80	90	100	110	120
ML7050 R110M Head (110 Micron)	80	90	100	110	120	130
ML7050 R130M Head (130 Micron)	100	110	120	130	140	150

(Based on average cornea thickness)

- Using a Calibrated LASIK Blade (CLB®) is the only way to achieve an accurate flap thickness with the ML7050 M Heads.
- MED-LOGICS also offers a 24 hour turn around repair service for the M2 handpiece.