


PRPS

REF	mini bags per kit	mini bag filling volume	pieces per box
PRPS 5	5	0 - 10 ml	25
PRPS 10	10	0 - 10 ml	25



PRPS 1

REF	mini bags per envelope	mini bag filling volume	pieces per box
PRPS 1	3	0 - 10 ml	25



PRPS 1C

REF	micro bags per envelope	micro bag filling volume	pieces per box
PRPS 1C	3	0 - 1 ml	25



KIT FOR
THE PREPARATION OF
ALIQUOTS OF PRP
AND PLATELET GEL
WITH FRESH
HEMOCOMPONENTS

PRP SECURE

BIOMED DEVICE S.R.L.
via Ippolito Nievo, 89
41124 Modena - Italia

telefono +39.059.343929 - fax +39.059.4550010
e.mail: info@biomeddevice.it
www.biomeddevice.it



PRACTICAL QUICK SAFE PRP SECURE

KIT FOR THE PREPARATION OF ALIQUOTS OF PRP AND PLATELET GEL


WITH FRESH AUTOLOGOUS HEMOCOMPONENTS
SYSTEM FOR HEMOCOMPONENTS PREPARATION
AND PRESERVATION.

PRPS is the results of the experience and long-term research activities of Biomed Device PRPS allows for the preparation of aliquots of fresh hemocomponents (PRP, CP, Platelet Gel).

It is a closed-circuit, pre-connected and sterile that offer maximum security and it is very practical to use. PRPS can be used in unprotected environments area also. **Every technical detail has been especially designed to fulfil the requirements of hospital workers in the most simple, fast and safe manner.**

DO YOU NEED SMALL DOSES?

Choose the single configuration PRPS 1 or PRPS 1C.



product
in envelopes
containing
3 mini
bags

lot complete
traceability

freeze
until
-80°C

injection
port



CHARACTERISTICS

- PRPS: ramp configuration, from 5 or 10 mini bags with injection port for activation or withdrawal. Filling volume for mini bag up to 10 ml.
- PRPS 1: single 10 ml mini bag.
- PRPS 1C: single 1 ml micro bag.
- Filling line with luer lock and valve for transfer from the bag. Filling line features diameter and materials that are compatible with the most commonly used sealing/sterile connection systems.
- Materials: EVA - PVC DOP / DEHP FREE - PE - PP - PC - ABS - PTFE. Latex-free device.
- Registered in the Directory of Medical Devices of the Italian Ministry of Public Health. CE 0051 Certification.
- Patent pending: all PRPS line products are protected by international industrial patents.

