

purple surgical

Operating Theatre Essentials

Incorporating Safe Sharps Management



FROM THE UK'S LEADING MANUFACTURER OF LAPAROSCOPIC INSTRUMENTS AND DEVICES

Sharps Safety Station™

The unique, patented Purple Surgical Sharps Safety Station™ range provides the ultimate in sterile sharps organisation, counting and disposal of contaminated sutures, blades and needles. This extensive range incorporates a number of unique safety features, which significantly reduce the risk of inadvertent sharps injuries.

Features and benefits include:

- + Patented, quick, simple and safe universal blade remover
- + Unique, convenient integral scalpel shield for safe and simple scalpel storage and retrieval
- + Unique locking mechanism prevents inadvertent opening whilst enabling re-counts if required
- + Unique, secure, three hinge design prevents box opening during transportation
- + Deep overlapping edges prevent accidental sharps migration
- + Four size options to meet all surgical requirements
- + Wide choice of formats, including foam blocks, magnetic surfaces, adhesive surfaces and foam strips, to meet customer preference and clinical indication
- + International bio-hazardous waste colour identification




Product	Code	Description	Size	Pack Qty
Sharps Safety Station™	PS2535	Double Magnet	Medium	100 / Box
Sharps Safety Station™	PS2560	Double Magnet	Large	100 / Box
Sharps Safety Station™	PS2635	Magnet and Foam Strips	Medium	100 / Box
Sharps Safety Station™	PS2615	Magnet and Foam Strips	Large	100 / Box
Sharps Safety Station™	PS2925	Magnet and Foam Block	Medium	100 / Box
Sharps Safety Station™	PS2940	Magnet and Foam Block	Large	25 / Box
Sharps Safety Station™	PS2922	Adhesive and Foam Block	Medium	100 / Box

Sharps Count

The Purple Surgical Sharps Count is a rigid plastic, puncture-proof box for the organisation, counting and disposal of contaminated sutures, blades and needles.

Features and benefits include:


- + Strong rigid construction prevents accidental penetration by contaminated sharps
- + A clear lid option allows convenient sharps visualisation and counting
- + Unique, locking mechanism prevents inadvertent opening whilst enabling re-counts if required
- + Unique, secure, three hinge design prevents box opening during transportation
- + Deep overlapping edges prevent accidental sharps migration
- + International bio-hazardous waste colour identification



Product	Code	Description	Size	Pack Qty
Sharps Count	PS3100	Magnet with Clear Top	Medium	100 / Box
Sharps Count	PS3110	Magnet with Clear Top	Large	100 / Box
Sharps Count	PS2910	Magnet and Foam Block	Small	100 / Box

Puncture-Proof-Pad

The Purple Surgical Puncture-Proof-Pad is a unique rigid plastic, puncture-proof-pad for the simple, safe organisation, counting and disposal of contaminated sutures, blades and needles.



Features and benefits include:

- + Strong rigid construction prevents accidental penetration by contaminated sharps
- + Wide adhesive strip ensures secure sharps retention
- + Secure locking mechanism prevents inadvertent opening
- + Overlapping edges prevent sharps migration
- + An optional 'fold flat' edge allows the safe accommodation of hypodermic needles and scalpel handles when required
- + International bio-hazardous waste colour identification

Product	Code	Description	Pack Qty
Puncture-Proof-Pad	PS2002	Adhesive and Foam	50 / Box

Purple Turtle

The Purple Turtle is a unique cost effective puncture resistant pad for the simple, safe organisation, counting and disposal of contaminated sutures, blades and needles.

Features and benefits include:

- + Tough puncture resistant plastic offers excellent protection against the risk of accidental sharps penetration
- + Wide adhesive strip ensures secure sharps retention
- + Overlapping edges prevent sharps migration
- + A castellated edge allows the safe accommodation of hypodermic needles




Product	Code	Description	Pack Qty
Purple Turtle	PS2001	Adhesive and Foam	60 / Box

Sharps Pad

The Purple Surgical Sharps Pad is a simple, low cost cardboard backed foam pad which offers the minimum requirement for the organisation, counting and disposal of contaminated sutures, blades and needles.

Features and benefits include:

- + Purpose designed
- + Convenient
- + Low cost



Product	Code	Description	Pack Qty
Sharps Pad	PS2000	Adhesive and Foam	200 / Box

Scalpel Shield

The Purple Surgical Scalpel Shield is a purpose designed device for the safe, secure storage of up to three scalpels in the sterile field.

Features and benefits include:

- + Conveniently stores handles in the 'upright' position ensuring safe and simple placement and retrieval
- + Shrouded blade design prevents accidental personal injury and damage to drapes
- + Adhesive tab secures and stabilises device in the sterile field



Product	Code	Description	Pack Qty
Scalpel Shield	PS4593	Scalpel Shield	100 / Box

Mag-Mat

The Purple Surgical Mag-Mat is a reusable magnetic mat that provides a secure 'hands free' transfer zone for metallic instruments, thus reducing the risk of third party inflicted sharps injuries.

Features and benefits include:

- + Strong ceramic magnetic surface ensures secure retention of metallic instruments
- + Will not magnetise instruments
- + Autoclavable
- + High quality, reinforced mesh design ensures a cost effective, long-lasting device
- + Large (30cm x 40cm) surface allows for containment of multiple instruments
- + Waterproof construction minimises the risk of any pathogen cross infection




Product	Code	Description	Pack Qty
Mag-Mat	PS4116	Reusable Mag-Mat 30cm x 40cm (12" x 16")	1 / Each

Magnetic Drapes

Purple Surgical Magnetic Drapes are lightweight, single-use and provide a secure 'hands free' transfer zone for metallic instruments, thus reducing the risk of third party inflicted sharps injuries.

Features and benefits include:

- + Strong magnetic surface ensures secure retention of metallic instruments
- + Convenient single-use device
- + Lightweight construction makes the drape particularly suitable for young, elderly or frail patients
- + Waterproof construction minimises the risk of any pathogen cross infection
- + Available in two sizes




Product	Code	Description	Pack Qty
Magnetic Drape	PS4121	Single-use Magnetic Drape 25cm x 40cm (10" x 16")	25 / Box
Magnetic Drape	PS4122	Single-use Magnetic Drape 40cm x 50cm (16" x 20")	25 / Box

Light Handle Cover

The Purple Surgical Light Handle Cover is a single-use, sterile cover that fits easily and securely over our universal operating theatre light handle, which is supplied on a 'free loan' basis.

Features and benefits include:

- ✦ Sterile, single-use convenience
- ✦ Gives control of light direction to surgeon and/or scrub personnel
- ✦ Securely retained by the unique 'rigid' flange and ribs
- ✦ Full range of universal handle adaptors ensure compatibility with all theatre lights
- ✦ Ability to conveniently 're-cover' if contaminated
- ✦ Available in single, double or triple packs for maximum convenience and minimal waste



Product	Code	Description	Pack Qty
Light Handle Covers	PS4611	Single pack	100 / Box
Light Handle Covers	PS4612	Dual pack	100 / Box
Light Handle Covers	PS4613	Triple pack	100 / Box

Skin Marker Pens

The Purple Surgical Skin Marker Pens allow for the precise marking on skin prior to and during surgery.

Features and benefits include:

- + Sterile, single-use convenience
- + Blue ink is quick drying and does not run
- + Enclosed flexible ruler ensures accurate measurement, even over contoured surfaces
- + Large ink reservoir
- + Full range of surgical applications due to the availability of standard, fine and dual tip options



Product	Code	Description	Pack Qty
Skin Marker Pens	PS3150	Skin Marker Standard Tip with Flexible Ruler	25 / Box
Skin Marker Pens	PS3150L	Skin Marker Standard Tip with Flexible Ruler and Label	25 / Box
Skin Marker Pens	PS3151	Skin Marker Fine Tip with Flexible Ruler	25 / Box
Skin Marker Pens	PS3151L	Skin Marker Fine Tip with Flexible Ruler and Label	25 / Box
Skin Marker Pens	PS3158	Skin Marker Dual Tip with Flexible Ruler	25 / Box
Skin Marker Pens	PS3158L	Skin Marker Dual Tip with Flexible Ruler and Label	25 / Box

Swab Counting Bags

The Purple Surgical Swab Counting Bags allow for the safe, effective, convenient counting, handling and disposal of contaminated swabs.

Features and benefits include:

- + Safe, hygienic system minimises the risk of pathogen inhalation
- + Clear construction allows visualisation for rapid, accurate counting
- + Large or small swab counting options
- + Enables simple swab weight measurement for blood loss analysis
- + Courtesy opening tab minimises handler's exposure to contaminants
- + Available with standard or extended width swab racks



Product	Code	Description	Pack Qty
Swab Counting Bags	PS3050	Swab Bags, Non Sterile	250 / Box
Swab Rack	PS3055	Swab Rack, Reusable	1 / Each
Swab Rack Extension Bar	PS3056	Swab Rack Extension Bar, Reusable	1 / Each

Abrasive Tip Cleaner

The Purple Surgical Abrasive Tip Cleaner allows for the effective removal of charred debris from electrodes.

Features and benefits include:

- + Large 5cm x 5cm effective area
- + Extended release liner enables convenient removal with gloved hands
- + Radio-opaque marker ensures patient safety through rapid detection
- + Coated abrasive surface ensures no flaking of contaminated abrasive material during use



Product	Code	Description	Pack Qty
Abrasive Tip Cleaner	PS5204	5cm x 5cm	100 / Box

Electrosurgical Pencils

The Purple Surgical Electrosurgical Pencils are available as hand operated or foot operated, with a full range of accessories.

Features and benefits include:

- + Sterile, single use convenience
- + Comfortable, ergonomic design
- + Tactile cut and coag button action
- + Universal ESU fitting
- + Extensive range of electrode accessories



Product	Code	Description	Pack Qty
Electrosurgical Pencil	PS5555	Handswitch Operated Pencil	25 / Box
Electrosurgical Pencil Accessories	PS5003	Ball Electrode	50 / Box
Electrosurgical Pencil Accessories	PS5004	Extended Ball Electrode	50 / Box
Electrosurgical Pencil Accessories	PS5005	Flat Blade Electrode	50 / Box
Electrosurgical Pencil Accessories	PS5006	Extended Flat Blade Electrode	50 / Box
Electrosurgical Pencil Accessories	PS5007	Needle Electrode	50 / Box

Vessel Loops

Purple Surgical Vessel Loops allow for the atraumatic retraction, occlusion and identification of arteries, veins, nerves and ureters.

Features and benefits include:

- + Soft, strong, atraumatic silicone construction
- + Oval cross section ensures even dispersion of pressure
- + Radio-opaque material ensures patient safety through rapid location
- + Colour coded for use with arteries, veins, nerves and ureters
- + Available in two sizes – mini and maxi



Product	Code	Description	Pack Qty	Pack Qty
Silicone Vessel Loops	PS3001	Mini	○ White	10 / Box
Silicone Vessel Loops	PS3002	Mini	● Red	10 / Box
Silicone Vessel Loops	PS3003	Mini	● Blue	10 / Box
Silicone Vessel Loops	PS3004	Mini	● Yellow	10 / Box
Silicone Vessel Loops	PS3201	Maxi	○ White	10 / Box
Silicone Vessel Loops	PS3202	Maxi	● Red	10 / Box
Silicone Vessel Loops	PS3203	Maxi	● Blue	10 / Box
Silicone Vessel Loops	PS3204	Maxi	● Yellow	10 / Box

Golden Retriever

The Purple Surgical Golden Retriever is a reusable device for the quick, simple and safe retrieval of metallic sharps that have been inadvertently dropped onto the operating room floor.

Features and benefits include:

- + Strong magnetic rollers ensure safe, 'hands free' sharps retrieval
- + Robust design ensures a cost effective, long lasting device
- + Long handle enables retrieval without the need for bending



Product	Code	Description	Pack Qty
Golden Retriever	PS4000	Reusable Magnetic Retriever	1 / Each

Instrument Pouch

The Purple Surgical Instrument Pouch enables the secure, convenient housing of instruments close to the surgical site, thus minimising the necessity for hand-to-hand transfer.

Features and benefits include:

- + Adhesive strip securely attaches the pouch to the patient drape
- + Deep pockets ensure the safe retention of all standard length minimally invasive instruments
- + Instrument separation is assured by the three pocket design
- + The clear material design ensures simple, quick instrument identification and retrieval



Product	Code	Description	Pack Qty
Instrument Pouch	PS4050	Instrument Pouch	25 / Box

Patient Eye Shields

The Purple Surgical Patient Eye Shields effectively cover and protect eyes from dust and dirt

Features and benefits include:

- + Universal design fits both left and right eyes, thus minimising inventory
- + Transparent material allows limited vision if required
- + High quality, smooth, ergonomic design for optimum patient comfort



Product	Code	Description	Pack Qty
Patient Eye Shield	PS5666NS	Universal Patient Eye Shield - Non Sterile	300 / Bag

Purple Surgical is the leading independent British manufacturer of high quality, cost-effective single-use surgical instrumentation and devices, partnering with global healthcare providers for positive patient outcomes.

Contact Us

Tel: +44 (0)1923 839333

Email: mail@purplesurgical.com

www.purplesurgical.com



Great British Manufacturing

Purple Surgical UK Ltd

Purple Surgical International Ltd

2 Chestnut House, Farm Close, Shenley, Herts. WD7 9AD, England